

# Solutions@ Mecmesin

## Lighting Weld Test

### Specification

A simple low-cost fixture was specified by the customer for the quality control laboratory. This was to hold a lighting assembly in a test where the welded contacts were pulled off.

A motorised test system was specified as the test needed to be repeatable and accurate.

### Solution

A multi-purpose small wedge grip was used to hold both electrical contacts. The bulb assembly was held inside a mounting fixture attached to the motorised test stand.

To prevent point contact and possible damage to the bulb assembly, 2 pegs with rubber sleeves were used to secure it. The contacts were pulled to destruction at a test speed of 50mm/min with peak load being captured at fast rate of 1200Hz to ensure a high degree of repeatability between tests. A cable between the AFG and VersaTest stand was used to automatically reverse the crosshead of the stand to its start position after the contacts were broken



### System

- VersaTest Motorised Test Stand
- 1000N AFG
- Bulb Assembly Mounting Fixture
- Small Wedge Grip

### Supplied to

GE Thorn Lighting, UK

Mecmesin Limited  
Newton House, Spring Copse Business Park, Slinfold,  
West Sussex, United Kingdom, RH13 0SZ

e: [sales@mecmesin.com](mailto:sales@mecmesin.com) t: +44 (0) 1403 799979 f: +44 (0) 1403 799975 www.mecmesin.com

## electrical & electronics industry

